

24/7 EEG™
SubQ

Specifications & Dimensions



EEG specifications

Sampling frequency	207 Hz
Bandwidth	0.5-48 Hz
Dynamic range	± 350 μV
Noise level	<1.75 μV_{rms}
Memory	+30 days
Data format	EDF+





Core features

- 2 channel EEG recordings
- 1 month data storage capacity
- 24 hour battery life

Direct user feedback, via

- Light
- Sound
- Vibration

Attachment Magnet





Core features






24/7 EEG™
SubQ

For more info
please contact

UNEEG™ medical A/S
UNEEG@uneeg.com





247™

EEG SubQ